

1. Restore natural knee function


Empowering Kinematic Restoration of the
Knee to restore natural Knee function

The Challenge

Help you to do better for your patients.

At 6 and 12 months postoperative, ~30% to 60% of patients report limitations in meeting functional expectations, like getting on/off a toilet or walking up/down stairs^{1*}

Pre-op Expectations Met/Exceeds at 6 & 12 months Post-op¹




* Study conducted in USA

** Study conducted in Canada

J&J MedTech

80% of TKA patients do not experience a natural joint vs. 42% of THA patients^{51**}

Patient Perception of TKA and THA²



Orthopaedics

Kinematic Restoration

Kinematic Restoration is the art and science of restoring knee function.

- From pre-op through post-op, from incision to closure—it provides a **holistic approach** that works in harmony with kinematic and inverse kinematic alignment philosophies
- The VELYS™ Robotic-Assisted Solution equips surgeons with **real-time** insights to optimize implant placement, preserve soft tissue envelope, predict joint stability and work toward returning knee function.^{5-7,42}
- ATTUNE™ Knee is verified to work with kinematic alignment techniques within tested boundaries.^{3-6*}

The powerful combination can reduce soft tissue releases to improve function and may accelerate patient recovery.⁷


*Based on in vitro studies.



Johnson & Johnson MedTech Knee Platform Solution

Kinematic Restoration


Surgeon driven
landmark acquisition




Surgeons control the 3D
surgical plan (e.g. rKA, iKA, MA)



Pre-resection soft tissue
balance check




Surgeon controlled,
Robotic-Assisted
resection



Immediate accuracy
verification & ability to
adjust the plan



Implantation of
Kinematically-advanced
ATTUNE Knee



Stability

Multiple studies have demonstrated **ATTUNE to exhibit smooth roll-back without abrupt change at mid-flexion**⁸⁻¹¹



Satisfaction

ATTUNE Knee is delivering better patient satisfaction than the class of other knees.^{12¶}

90%


Of patients are feeling better following TKA with ATTUNE Knee.^{1§}

88%

Of patients are satisfied with the results of their ATTUNE Knee TKA surgery.^{2¶}

Improved Accuracy

Studies demonstrate VELYS Robotic-Assisted Solution has improved implant alignment accuracy compared to Manual.^{13,14 †}



Fewer Soft-Tissue Releases

The VELYS Robotic-Assisted Solution showed significantly lower incidence of intentional soft tissue release compared to manual TKA.⁷



Improved Patient Function & Outcomes

After one year of follow-up, patients treated with the VELYS™ Robotic-Assisted Solution showed statistically significant improvements ($p<0.05$) in their Oxford Knee Score, WOMAC, Forgotten Knee Score, KOOS ADL, Activity Pain, and Range of Motion compared to those treated with the BRAINLAB KNEE3 Navigation System^{10*}



P<0.05
Oxford, WOMAC,
KOOS, FKS
vs. KNEE3 Navigation






¶According to the NJR ISR for ATTUNE Cemented, 88% of patients reported satisfaction in responding to the question "How would you describe the results of your operation?"
†Product may not be commercially available in all markets. This content may not be used externally in those markets where regulatory approval has not been granted for all the products referenced.
*In a prospective study of 90 TKA patients.
§ According to a Patient Satisfaction PROMS analysis 90% reported improvement in responding to the question "Overall, how are your problems now, compared to before your operation".

Limitations

How can you do better for your patients?




Implant

1. Research has shown that legacy **Designs can contribute to mid-flexion Instability** which can lead to a lack of patient confidence in the knee.¹⁶
2. The natural knee is seen to be very stable medially and roll back laterally, **medial instability has been associated with worse outcomes.**¹⁷
3. **Patella Femoral tracking can be compromised** when performing patient specific techniques.¹⁸





Technique

1. Mechanical alignment is the current gold standard tech, but **only 0.1% of patients have neutral femurs and tibias.**¹⁹
2. Mechanical alignment **requires soft tissue releases to balance the knee** in a high proportion of cases.²⁰
3. Implants designed and **tested under mechanical alignment conditions. Very limited published implant testing under alternative alignment conditions** which may influence wear and fixation performance.²¹




Technology

1. Manual instruments have **reduced accuracy** when implementing patient specific techniques **and lack predictive data on balancing.**²²
2. Navigation systems require **positioning and pinning of cutting blocks** adding additional steps to the procedure.²³
3. CT based Robotic systems are **reliant on technician defined landmark acquisition.**²⁴



J&J Knee Platform Solutions


Clinically empowering Kinematic Restoration of the Knee to restore natural Knee function



ATTUNE™ Knee System Implant



1. The ATTUNE GRADIUS™ Curve is the gradually reducing femoral radius designed to **provide a smooth transition from stability to mobility** through a patient's range of motion.⁸⁻¹¹
2. ATTUNE Medial Stabilized Knee System
 - Reduces A/P motion without abrupt direction changes
 - **Increases conformity on the medial compartment**
 - Lateral compartment is flatter to encourage natural femoral rollback^{8,11,25,26}
3. ATTUNE Knee Advanced Patellofemoral Articulation geometry **accommodates variations in trochlear orientation required for kinematic philosophies**^{27, 43-49}






J&J MedTech

Patient Specific Alignment Technique




1. Alternative alignment techniques such as rKA and iKA are capable of **restoring the alignment and joint line obliquity for a high proportion of the population**^{28,29}
2. Alternative alignment techniques such as rKA and iKA **reduce the need for soft tissue releases and adjustments**.²⁸⁻³⁰
3. ATTUNE has approved techniques supporting alternative alignment techniques and have published evidence on performance of ATTUNE in alternative alignments.³¹⁻³³



VELYS™ Robotic-Assisted Solution Technology




1. Parameters are easily adjusted through the VELYS Robotic-Assisted Solution planning screen, helping surgeons to **personalize alignment and balance relative to soft tissues**.^{32, 50*}
2. VELYS Robotic-Assisted Solution maintains the saw cut plane to help execute **precise, reproducible, surgeon-controlled cuts without cutting blocks**^{14*}
3. **Surgeons dictate** the selection of the critical points and define the 3D surgical plan³²



*Product(s) may not be commercially available in all markets. This content may not be used externally in those markets where regulatory approval has not been granted for all the products referenced.

ATTUNE Evidence


Clinically empowering Kinematic Restoration of the Knee to restore natural Knee function




ATTUNE™ Knee System Implant



- Multiple studies have demonstrated ATTUNE GRADIUS™ Curve to exhibit smooth roll-back without abrupt change at mid-flexion.⁸⁻¹¹
- Simulation has shown that the ATTUNE Medial Stabilized Knee System has lower overall medial condyle motion and does not see the abrupt changes in direction associated with implants based on a multi-radius sagittal geometry.²⁵
- A Study on the ATTUNE Knee with a patient specific technique found no influence of femoral rotation on PROMs and no patella complications, suggesting the patella femoral design is robust to positioning outside mechanical alignment.³⁴
- ATTUNE Knee is delivering better patient satisfaction than the class of other knees.^{12*}



Variable rotation of the femur does not affect outcome with patient specific alignment navigated balanced TKA
James Ferguson, PhD, The University of Texas at Austin



90%

Of patients are feeling better following TKA with ATTUNE Knee.¹¹

88%

Of patients are satisfied with the results of their ATTUNE Knee surgery.²⁷

∞ Shows the results of ATTUNE CR, ATTUNE MS, Persona MC, and Triathlon CS in both CR and CS applications. In CR application, medially, both the ATTUNE CR and MS implant constructs demonstrated minimal Femoral A-P translation, which stayed constant during most of the flexion. Both constructs experienced posterior femoral translation (rollback). The Persona MC experienced similar medial A-P translation as that of the ATTUNE CR and MS. However the A-P direction of translation reversed at about 25 degrees and 60 degrees of knee flexion. CS application, medially the ATTUNE MS construct experienced the least amount of total A-P translation, the ATTUNE MS, Persona MC and Triathlon CS (in that order). Conclusion: Demonstrates differences between the LP motion patterns, ligament force contributions and their strains for all simulated implant constructs in CR and CS application. The CS application, without PCL, the femur could slide anteriorly with flexion if the insert does not provide desired constraint to the anterior translation. Results: Demonstrated the effect of the size mis-matched Femurs in the Persona MC constructs. Trends indicated that the Femur-insert articulation for the smaller Femur size may offer less medial A-P stability that is achieved the larger sized Femur. The subsets of the results were normalized to their starting point at full extension and compared within designs. *According to the NJR ISR for ATTUNE Cemented, 88% of patients reported satisfaction in responding to the question "How would you describe the results of your operation?" *According to a Patient Satisfaction PROMS analysis 90% reported improvement in responding to the question "Overall, how are your problems now, compared to before your operation"

Technique





- iKA utilizing VELYS Robotic-Assisted Solution allows accurate reconstruction of joint line height & obliquity compared to MA.³⁵
- Avoiding soft tissue releases when implanting the ATTUNE™ Knee System has been associated with improved functional outcomes.⁷
- Published in vitro testing on the ATTUNE™ Knee Systems showed no increase in wear and equivalent cement fixation under simulated alternative alignment conditions as compared to mechanical alignment.^{33,36**}
- The use of iKA and KA have been associated with either equivalent or improved short term clinical and functional outcomes compared to mechanical alignment.^{37,38}

Table 2: Changes in joint line obliquity and joint line height

	Pre-op value	Post-op value	P-value
MA Group (n=101)			
Mean JLO, Degrees (SD)	2.94 (1.9)	2.31 (1.5)	0.004
Mean JLH, mm (SD)	40.6 (4.5)	40.6 (4.8)	0.89
iKA Group (n=101)			
Mean JLO, Degrees (SD)	2.43 (1.8)	2.30 (1.4)	0.57
Mean JLH, mm (SD)	41.2 (6.6)	42.4 (6.2)	0.17

JLO, Mechanical Alignment; iKA, Inverse Kinematic Alignment; JLH, Joint Line Obliquity; JLH, Joint Line Height; SD, Standard Deviation







**In vitro testing of components positioned with a simulated HKA of 3° Varus and a tibial joint line of 5° Varus

VELYS™
Robotic-Assisted Solution

Technology



- Studies demonstrate the use of VELYS Robotic-Assisted Solution results in improved implant alignment accuracy compared to manual.^{13,14 §}
- The VELYS Robotic-Assisted Solution showed significantly lower incidence of intentional soft tissue release compared to manual TKA.⁷
- VELYS Robotic-Assisted solution has demonstrated reduced navigation time and improved accuracy compared to Navigation.³⁹⁻⁴¹
- After one year of follow-up, patients treated with the VELYS™ Robotic-Assisted Solution showed statistically significant improvements (p<0.05) in their Oxford Knee Score, WOMAC, Forgotten Knee Score, KOOS ADL, Activity Pain, and Range of Motion compared to those treated with the BRAINLAB KNEE3 Navigation System^{15¶}



¶In a prospective study of 90 TKA patients.

§ Product may not be commercially available in all markets. This content may not be used externally in those markets where regulatory approval has not been granted for all the products referenced.

References

1. Tippet S. R. et al: Collecting Data with Palm Technology: Comparing Preoperative Expectations and Postoperative Satisfaction in Patients Undergoing Total Knee Arthroplasty, JBJS: December 1, 2010 - Volume 92 - Issue Supplement 2 - p 59-67 doi: 10.2106/JBJS.J.00827 Scientific Articles. (M_EM_DPS_JRKN_217451.1)
2. National Joint Registry for England, Wales, Northern Ireland and the Isle of Man. Implant Summary Report for DePuy ATTUNE Cemented. Web Site. www.ATTUNEvidence.com . Accessed 26 February 2025 (M_EM_DPS_JRKN_217461.1)
3. Johnson and Johnson MedTech. Vivo Knee Simulator Wear Test – Anatomical Alignment Attune TKA. 2017. Windchill #103384220. (M_US_DPS_JRKN_390282)
4. Johnson and Johnson MedTech. Attune Fixed Bearing: Intraoperative ROM Simulation and Pull-Off Test with Lipid and Marrow. 2017. Windchill #103397037.
5. Johnson and Johnson MedTech. Attune porous tibia micromotion with anatomic and mechanical alignment. 2021. Windchill #103839160. (M_US_DPS_JRKN_390282)(C_EM_DPS_JRKN_390281.EN)
6. Johnson and Johnson MedTech. Anatomic Alignment Surgical Technique and Patellar Tracking Validation Report. 2017. Windchill #103335944. (C_EM_DPS_JRKN_390281.EN) (M_US_DPS_JRKN_390282)
7. Spitzer A, Gorab R, Barrett W, et al. Robotic-assisted total knee arthroplasty reduces soft-tissue releases which improves functional outcomes: A retrospective study. Knee. 2024;49:52-61.(C_EM_DPS_JRKN_390281.EN)(M_US_DPS_DGSR_328959.1)
8. List E, Schutz R, Angst M, Ellenberger L, Datwyler K, Ferguson SJ. Videofluoroscopic Evaluation of the Influence of a Gradually Reducing Femoral Radius on Joint Kinematics During Daily Activities in Total Knee Arthroplasty, Journal of Arthroplasty, 2020; 35 (10), 3010-3030 (M_EM_DPS_JRKN_217508.1) (C_EM_DPS_JRKN_311748.EN)
9. Pfizner T, Moewis P, Stein P, et al. Modifications of femoral component design in multi-radius total knee arthroplasty lead to higher lateral posterior femoro-tibial translation. Knee Surg Sports Traumatol Arthrosc. 2018;26(6):1645-1655. doi:10.1007/s00167-017-4622-7 (M_EM_DPS_JRKN_217508.1)
10. Khasian M, Sharma A, Fehring TK, Griffin WL, Mason JB, Komistek RD. Kinematic Performance of Gradually Variable Radius Posterior-Stabilized Primary TKA During Various Activities: An In Vivo Study Using Fluoroscopy. J Arthroplasty. 2020;35(4):1101-1108. doi:10.1016/j.arth.2019.10.053 (M_EM_DPS_JRKN_217508.1)
11. Takagi H, Asai S, Sato A, Maekawa M, Kawashima H, Kanzaki K. Case series report of navigation-based in vivo knee kinematics in total knee arthroplasty with a gradually reducing femoral radius design. Ann Med Surg (Lond). 2017;17:33-37. Published 2017 Mar 29. doi:10.1016/j.amsu.2017.03.032 (M_EM_DPS_JRKN_217508.1)(C_EM_DPS_JRKN_311748.EN)
12. National Joint Registry for England, Wales, Northern Ireland and the Isle of Man. Implant Summary Report for DePuy ATTUNE Cemented. Web Site. www.ATTUNEvidence.com . Accessed 26 February 2025 (M_EM_DPS_JRKN_217461.1)
13. Alton TB, Severson EP, Ford MC, Lesko J, Leslie IJ (2025) VELYS robotic-assisted total knee arthroplasty: Enhanced accuracy and comparable early outcomes versus manual instrumentation during adoption. J Exp Orthop 12:e70163 (M_US_DPS_DGSR_105103.1)
14. Doan, G.W., et al., Image-Free Robotic-Assisted Total Knee Arthroplasty Improves Implant Alignment Accuracy: A Cadaveric Study. J Arthroplasty, 2022. 37(4): p. 795-801.
15. Clatworthy M. VELYS Robotic-Assisted Solution (VELYS Robotic-Assisted Solution) Vs. BrainLab Knee3 TKA - Early Outcomes Study. NZOA & AOA Combined Annual Scientific Meeting; 2022; Christchurch (M_US_DPS_DGSR_328488)
16. Clary, C.W., et al., The influence of total knee arthroplasty geometry on mid-flexion stability: an experimental and finite element study. J Biomech, 2013. 46(7): p. 1351-7
17. Tsukiyama H, Kuriyama S, Kobayashi M, et al. Medial rather than lateral knee instability correlates with inferior patient satisfaction and knee function after total knee arthroplasty. Knee.2017;24(6): 1478-1484. doi:10.1016/j.knee.2017.09.004
18. Lustig, S., E. Servien, and C. Batailler, How to optimize patellar tracking in knee arthroplasty? Orthop Traumatol Surg Res, 2023. 109(1S): p. 103458
19. Almaawi AM, Hutt JRB, Masse V, Lavigne M, Vendittoli PA (2017) The Impact of Mechanical and Restricted Kinematic Alignment on Knee Anatomy in Total Knee Arthroplasty. J Arthroplasty 32:2133-2140
20. Oussedik S, Abdel MP, Victor J, Pagnano MW, Haddad FS (2020) Alignment in total knee arthroplasty. Bone Joint J 102-B:276-279
21. Srivastava A, Lee GY, Stekllov N, Colwell CW, Jr., Ezzet KA, D'Lima DD. Effect of tibial component varus on wear in total knee arthroplasty. 5 presented at The Knee; 2012, 2012
22. Koutp, A., et al., Accuracy of Conventional Instrumentation is Dependent on Alignment Philosophy Using the Identical Surgical Technique in Total Knee Arthroplasty. J Knee Surg, 2024. 37(1): p. 20-25
23. DePuy Synthes, PATIENT SPECIFIC ALIGNMENT: TIBIA FIRST SURGICAL TECHNIQUE A BALANCED APPROACH USING BRAINLAB KNEE3 NAVIGATION.2019
24. Stryker,MAKO TKA™ Surgical Guide. 2016
25. DePuy Synthes. Forward dynamics modeling analysis: multi-patient verification report. 2021.Windchill. SEA103845908 (M_US_DPS_JRKN_311748)
26. Mane AM, Komistek RD, Heldreth MA. The soft tissue harmony with the knee motion in Cruciate Retaining Arthroplasty. ORS-2021. (C_EM_DPS_JRKN_311748.EN)
27. Mane, A., et al., Does Femoral Flange Design Influence Patellafemoral Mechanics in Patient Specific Alignment in ISTA 2023 (C_EM_DPS_JRKN_314345.EN)
28. Almaawi AM, Hutt JRB, Masse V, Lavigne M, Vendittoli PA (2017) The Impact of Mechanical and Restricted Kinematic Alignment on Knee Anatomy in Total Knee Arthroplasty. J Arthroplasty 32:2133-2140
29. Winnock de Grave P, Luyckx T, Van Criekinge T, Muller JH, Olivier B, Van Eecke E, et al. (2023) Inverse kinematic alignment accommodates native coronal knee alignment better in comparison to adjusted mechanical alignment and restricted kinematic alignment. Knee Surg Sports Traumatol Arthrosc 31:3765-3774
30. Oussedik S, Abdel MP, Victor J, Pagnano MW, Haddad FS (2020) Alignment in total knee arthroplasty. Bone Joint J 102-B:276-279

References

31. DePuy Synthes, ATTUNE® KNEE SYSTEM FEMUR FIRST ANATOMIC ALIGNMENT SURGICAL TECHNIQUE. 2017
32. Johnson & Johnson and its affiliates. VELYS™ Robotic-Assisted Solution for Total Knee. Surgical Technique, Femur First Approach (November 2024). Windchill Report 103744585 Rev J
33. Maag C, Cracaoanu I, Langhorn J, Heldreth M (2021) Total knee replacement wear during simulated gait with mechanical and anatomic alignments. Proc Inst Mech Eng H 235:515-52 (M_EM_OR_T_JRKN_396989)
34. Murgier J, Clatworthy M. Variable rotation of the femur does not affect outcome with patient specific alignment navigated balanced TKA. Knee Surg Sports Traumatol Arthrosc. 2022;30(2):517-526. <https://pubmed.ncbi.nlm.nih.gov/32783080/>. Accessed March 3, 2025. (M_EM_OR_T_JRKN_396987)
35. Shichman I, Hadad A, Brandstetter AS, Ashkenazi I, Warschawski Y, Gold A, et al. (2024) Better restoration of joint line obliquity in tibia first restricted kinematic alignment versus mechanical alignment TKA. Archives of Orthopaedic and Trauma Surgery;10.1007/s00402-024-05551-8 (M_EM_OR_T_DGSR_396986)
36. Maag C, Peckenpau E, Metcalfe A, Langhorn J, Heldreth M. INFLUENCE OF ALIGNMENT-SPECIFIC PHYSIOLOGICAL LOADING FOLLOWING INTRAOPERATIVE LIPID/MARROW INFILTRATION AND INTRAOPERATIVE MOTIONS UPON TIBIAL IMPLANT FIXATION. Orthopaedic Proceedings 2020;102-B(SUPP_1):32 (M_EM_OR_T_JRKN_396989)
37. Elbuluk AM, Jerabek S, Suhardi J, Sculco P, Ast M, Vigdorchik JM (2022) Head-to-Head Comparison of Kinematic Alignment versus Mechanical Alignment for Total Knee Arthroplasty. J Arthroplasty;10.1016/j.arth.2022.01.052
38. Winnock de Grave P, Kellens J, Tampere T, Vermue H, Luyckx T, Claeys K (2022) Clinical outcomes in TKA are enhanced by both robotic assistance and patient specific alignment: a comparative trial in 120 patients. Arch Orthop Trauma Surg;10.1007/s00402-022-04636-6
39. Ho JPY, Jagota I, Twigg J, Liu DW (2024) Robotic-Assisted Total Knee Arthroplasty Results in Shorter Navigation Working Time With Similar Clinical Outcomes Compared to Computer-Navigated Total Knee Arthroplasty. J Arthroplasty;10.1016/j.arth.2024.09.026 (M_EM_OR_T_DGSR_396988)
40. Bourgeault-Gagnon Y, Salmon L.J, Lyons MC (2024) Robotic-Assisted Total Knee Arthroplasty Improves Accuracy and Reproducibility of the Polyethylene Insert Thickness Compared to Manual Instrumentation or Navigation: A Retrospective Cohort Study. Arthroplasty Today 30 (M_EM_OR_T_DGSR_396988)
41. Doan GW, Van Avery A, Courtis RP, Leslie IJ, Hoeffel D, Clary C. Improvement in Primary Resection Accuracy with Image Free Robotic Assisted Total Knee Arthroplasty Compared to Computer Assisted Surgery (CAS). THE 21ST ANNUAL MEETING OF THE INTERNATIONAL SOCIETY FOR COMPUTER ASSISTED ORTHOPAEDIC SURGERY; 2022; Brest, France (M_EM_OR_T_DGSR_396988)
42. Huang P, Cross M, Gupta A et al. Early Clinical and Economic Outcomes for the VELYS Robotic-Assisted Solution Compared with Manual Instrumentation for Total Knee Arthroplasty. J Knee Surg 2024;37:864–872. (C_EM_DPS_JRKN_390281.EN)
43. Clary, C.W., Wright, A.P., Komosa, M.C., Maletsky, L.P. (2012b). The effect of patella medialization on patellofemoral kinematics after total knee replacement. 18th Congress of the European Society of Biomechanics (ESB); 29: 1262. (C_EM_DPS_JRKN_314345.EN)
44. Fitzpatrick CK, Clary CW, Rullkoetter PJ. The role of patient, surgical, and implant design variation in total knee replacement performance. Journal of Biomechanics. 2012;45:2092-2102. (C_EM_DPS_JRKN_314345.EN)
45. Barrett D, Wright, A., ATTUNE Knee System: GLIDERIGHT Articulation Whitepaper. (C_EM_DPS_JRKN_314345.EN)
46. Deacy J, Fitzpatrick C. Comparison of Natural and Unresurfaced Patellofemoral Mechanics. ORS 2012. Poster #1062. (C_EM_DPS_JRKN_314345.EN)
47. Clary CW, Fitzpatrick CK, Maletsky LP, Rullkoetter PJ, Probabilistic Modeling of Fem Component Overhang. ORS Poster #1988. 2012. (C_EM_DPS_JRKN_314345.EN)
48. Fitzpatrick C, Baldwin M. Identifying Alignment Parameters Affecting Implanted Patellofemoral Mechanics. J. of Orthopedic Research. JULY 2012. 1167-75. (C_EM_DPS_JRKN_314345.EN)
49. Johnson & Johnson MedTech (DePuy Synthes) IFU - ATTUNE Total Knee System - Fixed Bearing and Rotating Platform Bearing Knee. 0902-60-062_eIFU Rev. A. (C_EM_DPS_JRKN_314345.EN)
50. Johnson & Johnson and its affiliates. VELYS™ Robotic-Assisted Solution for Unicompartmental Knee Surgical Technique Guide. Rev B. May 27, 2025. Windchill #501965473.
51. Collins M, Lavigne M, Girard J, Vendittoli PA (2012) Joint perception after hip or knee replacement surgery. Orthop Traumatol Surg Res 98:275-280